Mathis
Labor-Heissluftofen
Typ LTH

Mathis
Mathis
Mathis
Mathis
Mathis

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This hot air oven is suitable for drying, hardening, baking and burn-in of applied primers, grounding agents, paint, lacquers or similar surface coating agents. The sample materials such as steel sheets, aluminium etc. can be coated according to the coil coating process by means of our reverse roll coater RRC-BW. This is a discontinuous unit which can duplicate coil coating applications of continuous production units. After the application the sample material is placed onto a holder and introduced into the dryer for the following thermal treatment.

A great number of important laboratory trials can be realised with this combination of application unit and drying oven. This avoids a more expensive and time consuming experiments with production ranges at a later stage.

When developing the concept of this dryer special care has been given to the absolute reproducibility of the parameters for the planned treatments. This is the only way to repeat the very same conditions of an earlier treatment.

The most important dryer components are:

Central control unit

Instead of a control box this apparatus has a compact control unit on which the following control parameters can be entered:

- temperature
- dwell time
- fan revolutions

The central control unit can be supplied with following options

- interface RS 232, for the online recording of the processing data
- autoprogrammer for 20 programs 20 steps each for the entry of controlled heating rates, holding times, fan speeds etc. This enables multiple treatments under various temperatures and dwell times.

Dryer box

By the use of a powerful fan in combination with a suitable, well designed air channel an excellent air and heat distribution could be achieved. The air channel is equipped with a regulation flap which guides the airflow either to the top or bottom only or to both sides at equal rates. The regulation flap can be set at any time from outside.

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Temperature measuring and heating


The temperature measuring is based on the resistance measuring system Pt 100. The temperature is controlled by an electronic measuring and regulating instrument with digital display. The heating elements are located above the air channel at the entrance of the fan. By means of an additional temperature sensor with corresponding regulator the maximum possible ambient temperature is limited in order to protect the heating elements.

The unit can be supplied with a non contact infrared radiation pyrometer for the measuring of the actual surface temperature of the sample material.

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Dwell time

This can be adjusted at the central control unit in minutes and seconds. The adjusted set value as well as the running down dwell time are displayed.

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Small space requirements

By eliminating the voluminous control box the space requirements could considerably be reduced. Nevertheless the electric components are well accessibly placed under the dryer chamber. The accessibility is guaranteed by easily removable side coverings.

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Transport frame

The entry of the material holder is activated by pressing a corresponding button and after running down of the selected dwell time the sample is delivered automatically.

Reliability

Several supervising systems are installed into this unit in order to increase the reliability, for instance the above mentioned temperature control for the heating. Furthermore the fan is equipped with a rotation control which cuts off the heating in case of problems with the fan motor or breakage of the fan belt.

Accessibility to the inside of dryer

The accessibility to the inside is guaranteed by removing the dryer cover. After that the whole inside is free for any check and maintenance.

Exhaust fan

An exhaust fan is available in case of heat treatments which cause considerable fumes and therefore require an exhaust fan.

The following accessories are available in order to adapt the hot air oven to your specific use:

Central control unit

- interface RS 232, for the online recording of the processing data
- auto programmer for 20 programs 20 steps each for the entry of controlled heating rates, holding times, fan speeds etc. This enables multiple treatments under various temperatures and dwell times.
- Programmer with weekly timer for
 - 4 ON/OFF switch cycles per day
 - 1 ON/OFF hand switch cycle (up to 5 days in advance)

Hot air oven

- Non contact infrared radiation pyrometer for the measuring of the actual surface temperature of the sample material
- Temperature increase up to 450°C
- Sieve frame sample holder with standing out cup screws
- Sieve frame sample holder with standing out cup screws as well as tilting device plus water trough for the fast cooling down of the sample after leaving the oven
- Lifting motor for the opening of the rear cover part in order to have a fast cooling down

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
- Exhaust fan, capacity max. 30 m³/h
- Foot switch for entry and delivery of the sample without hand operation of the key board
- Under carriage with rolls and brake

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LTH with tilting device plus water trough for the fast cooling down of the sample after leaving the oven

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Our manufacturing program furthermore includes:

- LABCOATER LTE-S, combination of drying chamber and coating device for the application onto textile materials, foils, paper etc.
- Roller application unit DWE for applications by means of engraved roller
- continuous laboratory and pilot coating ranges for smaller production widths up to approx. 60 cm
- laboratory apparatus for dyeing and finishing textiles and similar materials

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Werner Mathis AG
Werner Mathis U.S.A. Inc.
Concord N.C. 28026, U.S.A.